

Uracil DNA Glycosylase (UNG), heat-labile

Code No. 2820 **Size:** **200 U**
Conc.: **2 U/ μ l**

Description:

Uracil DNA Glycosylase (UNG) hydrolyzes N-glycosylic bonds between the deoxyribose sugars and the uracil bases in uracil-containing DNA leaving apyrimidinic sites in the DNA. UNG excises uracil from both single- and double-stranded dU-containing DNA but not from RNA.

Storage Buffer:

20 mM	Tris-HCl, pH 8.0 (at 4°C)
100 mM	KCl
1 mM	DTT
0.1 mM	EDTA
0.5% (v/v)	NP-40
0.5% (v/v)	Tween 20
50% (v/v)	Glycerol

Storage: -20°C

Source:

Produced in *E. coli* strain expressing a recombinant *Xiphophorus maculatus* UNG mutant.

Unit definition:

One unit of UNG is defined as the amount of enzyme required to digest 1 μ g of uracil-containing dsDNA at 25°C in 30 min.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Purity:

Endonuclease activity is not detected after incubation of 1 μ g pBR322 DNA with 50 units of this enzyme for 2 hours at 37°C as judged from agarose gel electrophoresis pattern.

Inactivation by heat:

The enzyme is completely and irreversibly inactivated by heat incubating at 50°C for 10 min.

Usage:

[Standard protocol for PCR Carryover prevention]

1. Prepare the following PCR reaction mixture.

<i>TaKaRa Taq</i> [™] Hot Start Version (5 U/ μ l)	0.25 μ l
10X PCR Buffer (Mg ²⁺ plus)	5 μ l
dU plus dNTP Mixture (Cat. #4035)	4 μ l
25 mM MgCl ₂	1.5 μ l
UNG (2 U/ μ l)	0.5 μ l
Primer 1	10 - 50 pmol
Primer 2	10 - 50 pmol
Template	1 - 5 μ l
Sterile purified water	up to 50 μ l

2. Incubate for 10 min at 25°C.
3. Incubate for 2 min at 95°C to heat-inactivate UNG.
4. Start PCR cycling program.

Regarding detailed protocol, please refer to the product manual for *TaKaRa* PCR Carryover Prevention Kit (Cat. #6088).

- * Suitable amount of the enzyme is dependent on the application. It is commonly in the range between 0.1 - 1 U/50 μ l of reaction volume.
- * This enzyme works in PCR and RT-PCR buffers commonly used.

TaKaRa Taq is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

